

The Millikan Oil Drop Experiment

University Physics Experiment Center UESTC



The Millikan Oil Drop Experiment



- Beginning in 1907, Millikan (won the 1923 Nobel Prize in physics) proceeded with the measurement of electronic charge.
- The experimental equipment is simple and effective.
- The idea and method are elegant and concise.
- The macroscopic mechanics model is used to study the quantum properties of the microscopic world.
- The data are accurate and stable, which is praised as the model of experimental physics.

Experiment content

1. Aim
2. Principle
3. Experiment set-up
4. Process and cautions
5. Data and discussion
6. Exercise

1 Aims

- (1) Study on the **Design** of Millikan Oil Drop Experiment.
- (2) By measuring the movement of the charged oil droplets in the gravitational field and the electrostatic field, the **charge discontinuity** is verified and the **fundamental charge value**.
- (3) Through the adjustment of experimental instruments, the selection of oil droplets, tracking and measurement, as well as experimental data processing, training **scientific experimental attitude**.

2 principle


2.1 $U=0$, drop with uniform speedForce balance $f = mg$

For: $6\pi\eta av = \frac{4}{3}\pi a^3 \rho g$

as
$$a = \sqrt{\frac{9\eta v}{2\rho g}}$$

Correct the viscous coefficient
$$\eta' = \frac{\eta}{1 + b/ pa}$$
2.2 $U \neq 0$ 时, Oil droplets still

$$qU / d = mg$$



$$q = \frac{18\pi d}{U \sqrt{2\rho g}} \left[\frac{\eta l}{t \left(1 + \frac{b}{pa} \right)} \right]^{3/2}$$

2.3 discussion

2.3.1 Several known charge amount of the oil droplets q_i , for calculate e

(1) Solve the greatest common divisor

*difficult

*for measurement error, **no** greatest common divisor exists!

(2) (assuming initial value and constant correction)

(a) assuming and initial value e_0

(b) $n_i = [q_i / e_0]$

(c) $e_i = q_i / n_i$

(d) $\bar{e} = \sum_m e_i / m$

Request!

⊙ Oil droplets with a smaller amount of electricity

⊙ Sampling oil droplets enough

(3) How to assume the initial value of e_0 is more reasonable ?

2.3.2 choose oil

Principle: less charged, the measurement is more accurate.

(1) Uniform speed: 2cm between 10s–30s

bigger, faster the oil drops

Too big oil, more error for measurement

Too small oil, Brownian motion is evident

(2) Balance voltage between 100V–300V

The equilibrium voltage depends on the amount of charge

Little the charged quantity, more balance voltage, vice versa

(3) Conclusion: charged quantity should be 2–10.

(4) Drop time and voltage should not be the same, avoiding the same charged quantity.

(5) Advice: first set voltage for 200V, inject the oil, and wait for moment, choose the suitable oil, which have slow motion for experiment.

2.3.3 How to Evaluate Error

Processing:

(a) Assuming current accepted value $e_0 = e_{\text{accepted}}$


(b) $n_i = [q_i / e_0]$

(c) $e_i = q_i / n_i$

(d) $\bar{e} = \sum_n e_i / n$

(e) $E = \frac{e - e_{\text{accepted}}}{e_{\text{accepted}}} \times 100\%$

3 Equipment set-up



4 Process and cautions

- (1) Adjust the camera angle, for scale line **horizontal and vertical**.
- (2) Sprayer must keep the **nozzle up or obliquely above**.
- (3) Adjust the CCD focal length, to ensure **clearly oil droplets**.
- (4) Control movement of oil, to achieve oil measured **several times**.
- (5) Checking, any oil droplets charged **more than 10e**, should be removed, re-select the appropriate oil droplets
- (6) Be **patient** in the measurement process.

5 Data and process

Oil droplets	Balance voltage $V(\text{v})$		Fall time $t(\text{s})$	
	Measurement	Averaged	Measurement	Averaged
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				